



Heater

Press to turn on/off the heat function while massage continues.

Power

Press the POWER button to turn on the massager with heat. Auto timer will shut it off after 15 minutes. Press the Power button anytime to turn the unit off.

DR(Direction)

Press the DR button to change and maintain direction of the massage roller. When the unit is turned on, the massage roller will automatically change its rolling direction every 60 seconds.

Special features

- Special design Rolling massage system offers deep tissue Shiatsu massage.
- Flex handles let you easily target the right spot and apply pressure to the massage with your hands.
- Easy to operate. Equipped with Overheat Protection Device & programmed with auto timer to ensure safety.
- Come with a removable cotton cover for softer skin touch.
- Quiet operation, no noise. Durable leather-like cover.

Massage System



Silicone



Twisting



Springer

※ Please refer to the pictures below for a full body massage.



1 Shoulder massage



2 Neck massage



3 Waist massage



4 Thigh massage



5 Abdomen massage



6 Foot massage



7 Back massage



8 Calf massage



9 Hip massage



10 Arm massage



11 Ankle massage



12 Control the massage strength by adjusting the distance between armrests