

Introduction

Healing Q

This advanced massager provides an immediate soothing Shiatsu massage at the comfort of your own place, either at home or in the office, or when you're resting in the vehicle after a long drive.

힐링안마기
Q

A Perfect Gift for You & Your Loved Ones

Thank you for choosing Healing Q.

As a leading shiatsu massager brand in Korea, Healing Q is a professional manufacturer of electric massagers, specializing in massagers with armrests.

In 15 years of continuous technical innovation, Healing Q has developed 4 generations of massaging systems, with Twisting Massage System & Silicone Massage System to be its latest invention.

Healing Q is an award-winning brand trusted by oriental consumers, where Shiatsu massage was originated from.




Features



1 BUTTON

- Single button function
 - Short press the button
 - turn on the power
 - Short press the button when the product is in operation
 - turn on/off the heater
 - Long press the button when the product is in operation
 - turn on/off the power



- Silicone massage ball system
- Closer to body skin and softer massage experience
- Changing rolling direction each second automatically
- Easy-operated single-button mode and adjustable heating function
- High-quality PU materials
- Overheating protection device




Full body massage




Silicone massage ball

The 4th generation silicone massage ball system




A

Upgraded the original 4 massage balls to 12 massage balls
Specially designed for silicone massage ball system



B

Dual effect of massage feeling
More comfortable massage experience by using silicone materials



C

Close friend to body
Silicone density of big massage balls: 50 small massage balls: 65

【Silicone massage ball】 Patented in China and South Korea

Cautions and Notice

Cautions:

1. Do not use on an infant or invalid or unconscious person.
2. Do not use the product on body wounds or injuries.
3. Do not disassemble or modify the product on your own.
4. Do not drop the product or apply force on it as this may cause damage to the parts.
5. Do not use the products while sleeping.
6. For safety reason , do not use the product while driving.
7. Do not use the product continuously for more than 20 minutes.
8. The massager may cause uncomfortableness to patients with diabetes, acute illness, malignant tumor, high fever, skin disease, etc. Consult your doctor before using the product.
9. Keep close supervision when the product is used by, on or near Children.
10. Please do not wash the product or operate the product in moist environment, as it may cause electrical shock.
11. Please do not leave the product unattended when it is plugged in.
12. Please examine the covering before each use. And discard the product if any place shows any signs of deterioration.

Notice:

1. Confirm the voltage specified on the adaptor before using it.
2. Start the product after making sure the correct connection to transformer.
3. When using on neck or abdomen, always use it with wearing on clothes in order to protect your body skin.
4. Please remember to wipe the massager with a soft, damp cloth instead of abrasive cleaners.
5. Place the massager in a safe, cool and dry place when not in use.
6. Keep away from all solvents and harsh detergents.

Notice and Specification

Trouble Shooting :

1. When the massager is unable to turn the power on, that may be cause by poor insertion of power supply cord. Please ensure the power cord is properly plugged into the power socket.
2. When the massager functions are not working, that may be caused by overheating due to long period of using. Please be noted that the product will stop working automatically after continuous usage to prevent it from overheating. Please wait for more than one hour and the product could be used normally.

Warranty:

We warrant this product against defects in materials and /or workmanship under normal use for a period of ONE year from the date of purchase by the original purchaser

If the problem occurs under normal use and within the duration period of manufacturer ' s limited warranty, the products is repaired or replaced with a new one by the manufacturer at its own options free of charge to the customers.

Specification :

Item Name: The cushion massager
Rated Voltage: DC 12V 2A
Power Consumption: 24W
Cover Material: Poly Urethane
Motion: Shiatsu with Heating
Packaging dimension: 33.5*16*17.5 CM